

# Microsart® Research Mycoplasma

## Mycoplasma Detection Kit for qPCR

### Benefits

- Easy to use
- Highest flexibility
- Maximum reliability



### Product Information

Microsart® Research Mycoplasma enables a fast and robust detection of Mycoplasma DNA in cell culture supernatants most applicable in research and development. Carefully selected primer|probe combinations are highly specific for a region within the 16 S rRNA gene of at least 110 Mycoplasma species.

## Working Principle

2 µL of sample material, e.g. cell culture supernatant, can be added directly to the PCR reaction tube. For the detection of Mycoplasma DNA a TaqMan® real-time qPCR is used. Depending on the sample matrix the Sartorius spin-column based DNA preparation can be performed prior to PCR analysis to increase sensitivity or prevent inhibition. 200 µL sample volume can be used as starting material for DNA preparation if using the Microsart® AMP Extraction kit. 2 µL of isolated DNA extract are amplified in a qPCR cyclor and the evaluation can be performed with the standard cyclor software.

## Applications

The Microsart® Research Mycoplasma real-time PCR protocol is especially designed for fast and reliable screening of cell culture supernatants most applicable in research and development, e. g. biotech and bio-pharmaceutical research and development, university and governmental research groups. It is used for direct detection of *Mollicutes* (*Mycoplasma*, *Acholeplasma*, *Spiroplasma*) in cell culture, cell culture media components and derived biologicals.

## Fast Result

Microsart® Research Mycoplasma utilizes real-time PCR (qPCR). The kit can be performed with any type of real-time PCR cyclor able to detect the fluorescence dyes FAM™ and ROX™. The detection procedure can be performed within 3 hours.

## Easy Handling

The kit contains all essential components in a ready-to-use master mix.

- Screening with a small sample volume
- Cost saving (in case there is no EP 2.6.7 compliance required)

## TaqMan® Probes

The application of TaqMan® probes adds specificity to the PCR detection system. Highly specific results are already generated during the cycling process – no subsequent melting curve analysis is needed.

## Contamination Prevention

The kit contains dUTP instead of dTTP, so the option is available to degrade amplicons from previous analyses by using uracil-DNA glycosylase (UNG). Thus, the occurrence of false-positive results can be minimized. UNG is not included in the kit.

## Summary

For scientists and lab technicians who need to screen cell culture supernatants for Mycoplasma DNA, Sartorius offers the Microsart® RESEARCH Mycoplasma Detection Kit.

# Technical Specifications

Each kit contains reagents for 25 reactions. The expiry date of the unopened package is specified on the package label. The kit components are stored at +2 to +8°C. After opening and rehydration the components need to be stored below -18°C. The LOT specific Certificate of Analysis can be downloaded from the manufacturer's website ([www.minerva-biolabs.com](http://www.minerva-biolabs.com)).

Kit Component	25 Reactions
Order No.	SMB95-1005
Mycoplasma Mix	1 × lyophilized
Rehydration Buffer	1 × 1.0 mL
Positive Control	1 × lyophilized
PCR grade Water	1 × 1.0 mL

# Ordering Information

## Mycoplasma Kits

Description	Quantity	Order No.
Microsart® Research Mycoplasma	25	SMB95-1005

## Accessories

Description	Quantity	Order No.
Microsart® AMP Extraction	50 extractions	SMB95-2003

## Related Products


Description	Quantity	Order No.
Microsart® AMP Mycoplasma	25	SMB95-1001
Microsart® ATMP Mycoplasma	25	SMB95-1003

**Germany**

Sartorius Lab Instruments GmbH & Co. KG  
Otto-Brenner-Strasse 20  
37079 Goettingen  
Phone +49 551 308 0

**USA**

Sartorius Corporation  
565 Johnson Avenue  
Bohemia, NY 11716  
Phone +1 631 254 4249  
Toll-free +1 800 635 2906

 For further information, visit  
[www.sartorius.com](http://www.sartorius.com)